

Disposable Needles



Made of high quality stainless steel
 With large ID and high flow
 Semi-transparent needle-hub for the convenient observation of blood reflux
 Delicately-designed needle-point with features of sharpness, smooth puncture, less damage to the tissue, less painful feeling to the patient
 Size distinguished by the needle-hub & primes color for the clear recognition
 The flexible supply form: in bulk or pack, sterilized or unsterilized
 Production of all kinds of special needles required by our customers
 Packing types: in bulk or in blister


Spec

Gauge : 17G, 18G, 19G, 20G, 21G, 22G, 23G, 24G
 25G, 26G, 27G, 28G, 29G, 30G, 31G
 Length : ½", 1" ¼" 1½"


Cannula



ISO Specification



RW- Regular Wall
 TW- Thin Wall
 ETW- Extra Thin Wall
 UTW- Ultra Thin Wall



Cannula for Disposable Needle, Insulin syringe, Pen needle, Blood collection needle, lancet, etc.
 Material : SUS 304
 Can make various sizes

GAUGE	Hub color	Outside Diameter.(mm)	Inside Diameter.(mm)		
	ISO6009		ISO	RW	TW
32R-0.23		0.229~0.241	0.089~0.104	0.105~0.124	0.125~0.145 UTW 0.146~
31R-0.25		0.254~0.267	0.114~0.124	0.125~0.145	0.146~0.175 UTW 0.176~
30R-0.30	Yellow(122C)	0.298~0.320	0.133~0.164	0.165~0.189	0.190~0.239 UTW 0.240~
29R-0.33	Red (1785C)	0.324~0.351	0.133~0.189	0.190~0.239	0.240~0.264 UTW 0.265~
28R-0.36	Blue-green (320C)	0.349~0.370	0.133~0.189	0.190~	-
27R-0.40	Medium grey (429C)	0.400~0.420	0.184~0.240	0.241~	-
26R-0.45	Brown(478C)	0.440~0.470	0.232~0.291	0.292	-
25R-0.50	Orange(1505C)	0.500~0.530	0.232~0.291	0.292	-
24R-0.55	Medium purple (199C)	0.550~0.580	0.280~0.342	0.343	-
23T-0.60	Deep blue (299C)	0.600~0.673	0.317~0.369	0.370~0.459	-
22T-0.70	Black (447C)	0.698~0.730	0.390~0.439	0.440~0.521	0.460~
21T-0.80	Deep green (3405C)	0.800~0.830	0.490~0.546	0.547~0.609	0.522~
20T-0.90	Yellow(122C)	0.860~0.920	0.560~0.634	0.635~0.686	0.610~
19T-1.10	Cream(155C)	1.030~1.100	0.648~0.749	0.750~0.849	0.687~
18T-1.20	Pink (231C)	1.200~1.300	0.790~0.909	0.910~1.040	0.850~
17T-1.40	Red - violet (528C)	1.400~1.510	0.950~1.155	1.156~1.243	1.041~
16T-1.60	White	1.600~1.690	1.100~1.282	1.283~1.383	1.244~
15T-1.80	Blue - grey (431C)	1.750~1.900	1.300~1.459	1.460~1.559	1.390~
14T-2.10	Pale green	1.950~2.150	1.500~1.599	1.600~1.726	1.560~
13T-2.40	Purple	2.300~2.500	1.700~1.955	1.956~	1.727~
12T-2.70	Pale blue	2.650~2.850	1.950~2.234	2.235~	-
11T-3.00	Green-yellow	2.950~3.150	2.200~2.463	2.464~	-
10T-3.40	Olive brown	3.300~3.500	2.500~2.818	2.819~	-

Certificates



FDA



CE0120



ISO13485



KGMP

SHINA CORPORATION

Insulin Syringe & Safety Syringe, Pen needle Medical Device MFG.
 CE0120, FDA510K, KFDA GMP, ISO 13485 Certified

www.shinacor.co.kr
 since 1990



Lifetime
 Companion...

We are an alternative source for quality life disability due to diabetes and are always with our customer as your lifetime companion.

SHINA Co., Ltd.

#1005, 28, HWANGSAEUL-RO, 200 Beon-GIL, BUNDANG-GU,
 SEONGNAM-SI, GYEONGGI-DO, 13595, KOREA. (OWNER'S TOWER, SUNAE-DONG)
 Tel: 82-31-711-8180 / Fax: 82-31-711-8190
 email : shinacor@chol.com

Diabetes Care

Insulin Syringes



Description

The high transparency of the barrel permits the end user a strict control on the foreign substance and liquid flow.
 The inerasable graduation permits an easy reading, 40 unit and 100 unit are available.
 The medical grade silicone oil is used as the lubricant the purpose of getting a smooth plunger movement.
 Sterilized by EO gas, non-toxic, non-pyrogenic, single use only.
 Material for barrel and plunger: medical grade PP(poly-propylene)
 Material for gasket natural rubber/latex-free
 Fixed needle, no dead space
 Ultra fine needle
 Individual pack : Blister/Poly bag
 Volume:0.3ml/0.5ml/1ml

Spec

Gauge : 29G, 30G, 31G, 32G, 33G

Length : 12.7mm, 8mm, 6mm



Pen needles



Description

- Tri-bevel needle tip for easier insertion
- Medical grade silicone surface treatment for advanced injection comfort
- Easy flow of Insulin without any impurities
- Extensive range of needle lengths suitable for all types of patients
- Compatible with all major brand Insulin delivery pens
- Sterile, Non-toxic, Non-Pyrogenic


Spec

Gauge : 29G, 30G, 31G, 32G, 33G, 34G

Length : 12mm, 10mm, 8mm, 6mm, 5mm, 4mm

Smart pen needle

Smart design to reduce the risk of bending and breaking



Conventional pen needle



Safety Syringes



Insulin Safety Syringes

Insulin safety syringe with permanently attached needle are designed to protect patients and healthcare professionals and minimize medication waste.

0.5mL 29G X 1/2" (12.7mm)

1.0mL 29G x 1/2" (12.7mm)

Turberculin Safety Syringes

Tuberculin safety syringe with permanently attached needle are designed to protect patients and healthcare professionals and minimize medication waste.

1.0mL 28G X 1/2" (12.7mm)

1.0mL 27G X 5/8" (15.9mm)

1.0mL 27G X 1/2" (12.7mm)

1.0mL 25G X 5/8" (15.9mm)

1.0mL 25G X 1" (25.4mm)

Specialty Safety Syringes LDS(Low Dead Space)

1.5mL - 22 G x 1/2" specialty safety syringe with permanently attached needle is designed to protect patients and healthcare professionals and minimize medication waste.

1.5mL 22G X 1 1/2" (38.1mm)

Safety Disposable Syringes

3mL

Needle : Convertible Various Sizes

Specialty Syringes



LDS(Low Dead Space)

Specialty syringe with permanently attached needle is designed to minimize medication waste.

1.0mL 22G X 1 1/2"(38.1mm)

1.5mL 22G X 1 1/2" (38.1mm)

3.0mL 22G X 1 1/2" (38.1mm)

